

datasheet

# pco.edge 5.5 CLHS

cooled sCMOS camera

resolution  
**5.5 MPixel**

pixel size  
**6.5  $\mu\text{m}$  x 6.5  $\mu\text{m}$**

interface  
**CLHS FOL**



low readout noise  
1.0 e<sup>-</sup> (med)

high speed  
100 fps

high dynamic range  
30 000 : 1

high resolution  
2560 x 2160 pixels

shutter modes  
rolling & global shutter,  
global reset

## technical data

### image sensor

	slow scan	fast scan
<b>sensor technology</b>	scientific CMOS (sCMOS)	
<b>color type</b>	monochrome color (bayer pattern)	
<b>resolution (horizontal x vertical)</b>	2560 px x 2160 px	
<b>pixel size (horizontal x vertical)</b>	6.5 μm x 6.5 μm	
<b>sensor size (horizontal x vertical)</b>	16.6 mm x 14.0 mm	
<b>sensor diagonal</b>	21.8 mm	
<b>shutter mode</b>	rolling shutter (RS) with selectable readout direction global reset (GR) / global shutter (GS)	
<b>modulation transfer function (theoretical max.)</b>	76.9 lp/mm	
<b>peak quantum efficiency</b>	60 % @ 600 nm (monochrome)	
<b>spectral range</b>	300 nm - 1100 nm (monochrome)	
<b>dark current (typ.)</b>	0.6 e <sup>-</sup> /pixel/s @ +7 °C sensor temperature (RS/GR) 0.9 e <sup>-</sup> /pixel/s @ +7 °C sensor temperature (GS)	
<b>fullwell capacity</b>	30 000 e <sup>-</sup>	
<b>readout noise (typ.)<sup>1</sup></b>	1.4 e <sup>-</sup> rms (RS/GR) 1.0 e <sup>-</sup> med (RS/GR) / / /	1.5 e <sup>-</sup> rms (RS/GR) 1.1 e <sup>-</sup> med (RS/GR) 2.5 e <sup>-</sup> rms (GS) 2.2 e <sup>-</sup> med (GS)
<b>dynamic range (intra-scene)<sup>2</sup></b>	30 000 : 1 (90 dB) (RS/GR) / /	27 000 : 1 (89 dB) (RS/GR) 13 500 : 1 (83 dB) (GS)

<sup>1</sup> The readout noise values are given as median (med) and root mean square (rms) values, due to the different noise models which can be used for evaluation. All values are raw data without any filtering.

<sup>2</sup> The dynamic range value is calculated with the median value of the readout noise and rounded.

### frame rate table

#### vertical resolution reduction in fps

@ scan rate	slow scan	fast scan	
@ shutter mode	RS/GR	RS/GR	GS
<b>2560 x 2160</b>	33	100	50
<b>2560 x 1024</b>	70	212	105
<b>2560 x 512</b>	140	422	208
<b>2560 x 256</b>	279	838	409
<b>2560 x 128</b>	550	1651	789

#### typical resolutions in fps

@ scan rate	slow scan	fast scan	
@ shutter mode	RS/GR	RS/GR	GS
<b>1920 x 1080</b>	67	201	100
<b>1600 x 1200</b>	60	181	90
<b>1280 x 1024</b>	70	212	105
<b>640 x 480</b>	150	450	222
<b>320 x 240</b>	297	893	436



**camera**

	slow scan	fast scan
<b>max. frame rate @ full resolution</b>	33 fps (RS/GR) /	100 fps (RS/GR) 50 fps (GS)
<b>exposure time range</b>		500 $\mu$ s - 2 s (RS) 10 $\mu$ s - 2 s (GR) / 10 $\mu$ s - 100 ms (GS)
<b>dynamic range A/D<sup>1</sup></b>	16 bit	
<b>conversion factor<sup>2</sup></b>	0.46 e <sup>-</sup> /DN	
<b>pixel rate</b>	200 MPixel/s (RS/GR)	572 MPixel/s (RS/GR/GS)
<b>region of interest (ROI)</b>	horizontal: steps of 16 columns (min. 64) vertical: steps of 1 row (min. 16)	
<b>binning</b>	horizontal: x2, x4 (sum) vertical: x2, x4 (sum)	
<b>non-linearity</b>	< 0.6 %	
<b>dark signal non-uniformity (DSNU)</b>	< 0.3 e <sup>-</sup> rms (RS/GR) /	< 0.3 e <sup>-</sup> rms (RS/GR) < 3.9 e <sup>-</sup> rms (GS)
<b>photo response non-uniformity (PRNU)</b>	< 0.34 %	
<b>cooling temperature image sensor</b>	+7 °C stabilized (up to +27 °C ambient temperature)	
<b>cooling method</b>	forced air optional: liquid cooling	
<b>trigger input signals</b>	external exposure start, external exposure control, sequence trigger, acquire enable	
<b>status output signals</b>	exposure, busy, line	
<b>input / output signal connectors</b>	SMA	
<b>time stamp</b>	in image (1 $\mu$ s resolution)	
<b>data interface</b>	Camera Link HS FOL	

<sup>1</sup> The high dynamic signal is simultaneously converted at high and low gain by two 11 bit A/D converters and the two 11 bit values are sophisticatedly merged into one 16 bit value.

<sup>2</sup> According to EMVA1288, the conversion factor equals the inverse of the system gain and can be operational mode dependent.

**quantum efficiency**



### general


power supply	24 VDC ( $\pm 10\%$ )
power consumption	max. 32 W
weight	0.8 kg air-cooled 1.05 kg liquid-cooled
dimensions (height x width x length)	76 mm x 70 mm x 122.5 mm
operating temperature range	+10 °C to +40 °C
storage temperature range	-10 °C to +60 °C
humidity range (non-condensing)	10 % to 80 % (recommended < 65 %)
certifications	CE, FCC, UKCA

### optical interface

direct mounting distance	11.1 mm ( $\pm 10\%$ )
lens mounting	C-mount, F-mount
optional lens mounting	TFL-mount
optional lens remote control (only air-cooled variant)	EF-mount, EF-S-mount

Configure your optical setup with our **MachVis Lens Selector** online tool.

### dimensions



outlines of pco.edge 5.5 CLHS air-cooled (all dimensions given in mm)

## software

### Your first choice is pco.camware:

Our main camera control software enables control of most camera settings and facilitates image acquisition and storage.

You can customize it exactly to your needs using different layouts, styles and features.

### You prefer to use a different software:



Our cameras integrate with a range of third-party software applications.

In microscopy we offer dedicated support for µManager, while ensuring compatibility with other software maintained by their providers.

### You want to create your own application:

We feature a wide range of software development kits (SDK) for various programming languages, such as C++, Python, C#, LabVIEW, Matlab, and Java.

If you are looking for more general SDKs, we present pco.sdk and pco.recorder, our low-level SDKs with C interface.



Our software is available for Windows and Linux platforms.

Visit our **website** for detailed information, installation guidance, and Github projects.

## areas of application

3D metrology | biochip reading | brightfield microscopy | calcium imaging | digital pathology | flow cytometry | fluorescence microscopy | fluorescence recovery after photobleaching (FRAP) | Förster resonance energy transfer (FRET) | high-content screening | high-speed brightfield ratio imaging | high-throughput screening | hyperspectral imaging | image intensifier imaging | industrial quality inspection | lightsheet fluorescence microscopy (LSFM) | lucky astronomy | ophthalmology | particle tracking velocimetry (PTV) | photovoltaic inspection | single molecule localization microscopy (SMLM) – PALM, STORM, dSTORM, GSDIM | spinning disk confocal microscopy | structured illumination microscopy (SIM) | super-resolution microscopy | total internal reflection fluorescence microscopy (TIRF) | wafer inspection

### ordering information

<b>pco.edge 5.5 CLHS</b>	85108072201	camera system, 2560 x 2160 pixel, monochrome, global and rolling shutter, CLHS interface, air cooling
<b>pco.edge 5.5 WAT CLHS</b>	85108072205	camera system, 2560 x 2160 pixel, monochrome, global and rolling shutter, CLHS interface, liquid cooling
<b>pco.edge 5.5 C CLHS</b>	85108072233	camera system, 2560 x 2160 pixel, color, global and rolling shutter, CLHS interface, air cooling
<b>pco.edge 5.5 C WAT CLHS</b>	85108072245	camera system, 2560 x 2160 pixel, color, global and rolling shutter, CLHS interface, liquid cooling


**pco.**<sup>®</sup>

address: Excelitas PCO GmbH  
Donaupark 11  
93309 Kelheim, Germany

phone: (+49) 9441-2005-0  
(+1) 866-662-6653  
(+86) 0512-6763-4643

mail: [pco@excelitas.com](mailto:pco@excelitas.com)

web: [www.excelitas.com/pco](http://www.excelitas.com/pco)



[excelitas.com](http://excelitas.com)

  
**excelitas**<sup>®</sup>