

user manual

pc_o. firmware programmer for GigE interface



pc_o.

Target Audience: This camera is designed for use by technicians, engineers, and scientists.

In case of any questions or comments, please contact us at PCO.



telephone	+49 (0) 9441 2005 50
fax	+49 (0) 9441 2005 20
email	info@pco.de
postal address	PCO AG Donaupark 11 93309 Kelheim, Germany

The cover photo shows an exemplary PCO camera system.
The lens is sold separately.

Copyright © 2016 PCO AG (called PCO in the following text), Kelheim, Germany. All rights reserved. PCO assumes no responsibility for errors or omissions in these materials. These materials are provided as is without warranty of any kind, either expressed or implied, including but not limited to, the implied warranties of merchantability, fitness for a particular purpose, or non-infringement. PCO further does not warrant the accuracy or completeness of the information, text, graphics, links or other items contained within these materials. PCO shall not be liable for any special, indirect, incidental, or consequential damages, including without limitation, lost revenues or lost profits, which may result from the use of these materials. The information is subject to change without notice and does not represent a commitment on the part of PCO in the future. PCO hereby authorizes you to copy documents for non – commercial use within your organization only. In consideration of this authorization, you agree that any copy of these documents, which you make, shall retain all copyright and other proprietary notices contained herein. Each individual document published by PCO may contain other proprietary notices and copyright information relating to that individual document. Nothing contained herein shall be construed as conferring by implication or otherwise any license or right under any patent or trademark of PCO or any third party. Except as expressly provided, above nothing contained herein shall be construed as conferring any license or right under any PCO copyright. Note that any product, process, or technology in this document may be the subject of other intellectual property rights reserved by PCO, and may not be licensed hereunder.

Released May 2016 © PCO AG

By Thomas Hartmann

pco.firmware programmer for GigE interface V1.01 © PCO AG, Germany


FIRMWARE PROGRAMMER FOR GIGE INTERFACE

1.1 INTRODUCTION

This manual describes the firmware update procedure of the GigE interface card in a **pco.dimax S/HD/HS** or a **pco.1200-4000** camera system.


In order to update the firmware, please connect the camera to the computer using the USB interface and disconnect GigE. Unzip the firmware package to a folder of your choice and start the pco.programmer.

1.2 FIRMWARE UPDATE



The camera must be connected using the USB interface. Connect only one camera when updating the firmware. Switch on the camera and wait until you hear a beep.



First the selection dialog is shown. Please select **USB**, in case it is not selected.




Press **Connect**. Now the status information should display a camera with its serial number.

Click **OK** to proceed to the programming step.

Select **Load File** and open the **ppg-file: FWGigE_Euros201_Notos108.ppg**





The pco.programmer dialog shows the current firmware version and the version of the loaded (new) firmware.

Always select **both** devices GigE μP and GigE FPGA.

Updates are available for **GigE interface microcontroller** (GigE μP) and **FPGA** (GigE FPGA).


Press **program**, which will first create a backup of your current firmware. The camera will **beep** several times during the update process. The update will start after the backup process is finished.

Wait until the complete process is finished. This can take **a few minutes**.

NOTICE

DO NOT UNPLUG

DO NOT UNPLUG AND/OR SWITCH OFF ANY OF YOUR DEVICES, NEITHER THE PC, NOR THE CAMERA, DURING THE UPDATE.



After a successful firmware programming process you should see two green bars inside the programming dialog and the message **FW updated successfully** will pop-up.

Press **ok** and then **quit** to leave the pco.programmer.

NOTE

If you hear two beeps something went wrong. Please try to repeat the firmware update.

In order to connect with GigE please **unplug the USB cable** and **connect GigE**.

After a few seconds you should hear **one beep**. This indicates that the update has been finished successfully and the camera is ready for use.

ABOUT PCO



pco.

In 1987, PCO was founded with the objective to develop and to produce specialized, fast and sensitive video camera systems, mainly for scientific applications. Meanwhile the product range of PCO cameras covers digital camera systems with high dynamic range, high resolution, high speed and low noise, which are sold in the scientific and industrial market all over the world.

Currently PCO is one of the leading manufacturers of scientific cameras. Worldwide representatives, together with our own sales department and technical support assure that we keep in touch with our customers and their needs. The actual wide range of specialized camera systems is the result of technical challenge and product specific know-how. A design according to advanced techniques, a high standard of production and strict quality controls guarantee a reliable operation of the cameras. Our own developments in conjunction with an excellent contact to leading manufacturers of image sensors ensure our access to state-of-the-art CCD and CMOS technology for our cameras.

Since 2001, PCO is located in its own facility building in Kelheim at the shore of the beautiful and international river Danube. Here in the county Bavaria, which is well known for its excellent support and conditions for high technology companies, we share the benefits of the simple access to high performance products and services in the surrounding area.

Kelheim itself is a historical town, first documented in 866. The small city is founded at the confluence of the Danube and the Altmühl, which has been converted into the Rhine-Main-Danube bypass channel for water transport. Located in Danube valley, it is the heart of a beautiful river and forest covered lime plateau landscape. It's landmark, the Hall of Liberation, was built by Ludwig I. in 1863 on the Mount Michael and is visible from all over the city and valley. The beautiful Danube Gorge, which is protected as natural monument since 1840, is located between Kelheim and the famous abbey Weltenburg.

pco.