

X-Cite® XYLIS II Control Options

Control Options

Manual Control – speedDIAL	<p>ON/OFF Speed sensitive intensity control 1 favorite intensity setting UV OFF Left and right-handed control option</p>
USB	<p>ON/OFF Intensity control Compatible with X-Cite® 120PC, X-Cite exacte and X-Cite 120LED drivers</p>
TTL	<p>ON/OFF</p>
Foot Pedal (optional)	<p>ON/OFF All X-Cite XYLIS II models have input port, foot pedal sold separately</p>